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Telefon 05661.71.3400
Fax 05661.71.3702

www.sartorius-bbi-systems.com

CERTOMAT® BS-T Incubation Shaking Box

Operating Manual

Introductory Note

This Operating Manual for the Incubator Shaking Cabinet CERTOMAT® BS-T refers to the design and equipment features of this device at the date of print given with the Rev. No. in the bottom line. The information included in this Operating Manual has been carefully estimated and checked. However, no guarantee can be given with respect to completeness and correctness of this information. Sartorius Stedim Biotech GmbH reserves the right to modify the equipment of its products as well as the documentation without notice.

Due to the continuous development of the units it may be possible that different building components or their use are not described. There may be information which does not apply to your device. If you find missing, misleading or wrong information or if you should wish additional information about different building components please send us a note.

The CERTOMAT® BS-T is just an example of laboratory equipment products from the Sartorius Stedim Biotech GmbH. Our range of equipment includes shakers and incubator shakers, centrifuges, ultrapure water systems, air samplers and fermenters/bioreactors.

You will find further information on our homepage www.sartorius-stedim.com

If you are interested in more information about the complete product program or about specific equipment, please contact us.

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany
Telephone +49.551.308.3118
Telefax +49.551.308.3918
info@sartorius-stedim.com
www.sartorius-stedim.com

EG-Declaration of Conformity

With the Statement of Conformity, which is attached or which accompanies the equipment, Stedim Biotech GmbH confirms the conformance of the CERTOMAT® BS-T with the specified guidelines.

General Safety Notes

- Laboratory seniors or the persons responsible for the use of the equipment must inform the personnel about the dangers at the operation of the CERTOMAT® BS-T. Furthermore they have to guarantee that the relevant safety conditions are observed. The personnel must have the necessary qualification for the operation of the equipment and its handling.
- Please take special care regarding the operating information. Improper use of the unit is thus avoided and you have full claim for guarantee. Please contact us if you have questions.

Steps to be carried out with special care or special aspects or special issues to be considered are marked as this paragraph.

Important safety notes in this manual are marked with this symbol. Such instructions apply in addition to the other safety instructions for the field of application and for the working place.

- Do not use devices with the wrong potential version.
- Disregard of the safety instructions may for example result in damages at the equipment or can cause other material and eventually also person damages.

Guide through this manual

This Operating Manual is organized systematically. Depending on your knowledge of the device you can start from the beginning or go directly to the chapter of your interest.

- Section 1 gives an overview on design and function of the equipment of the shaking box.
- Section 2 describes the preparatory measures which are required for installation at the working place. Furthermore the installation and the connection of the unit are explained.
- Section 3 includes the operating information. If the equipment is installed ready for use or if you are used to the CERTOMAT® BS-T, respectively, you can go straight on to this section.
- Section 4 is a summary of all important technical data and includes ordering information of optional equipment and accessories for the shaker.

Contents

Introductory Note

EG Declaration of Conformity

General Safety Notes

Guide through this manual

1	Design and Function	4	4	Technical Data, Ordering Information	39
1.1	Design and Construction Features	4	4.1	Technical Data	39
1.2	Available Equipment and Accessories	4	4.1.1	Mechanical Design	39
1.2.1	Shaking Drives	4	4.1.2	Mains Supply	39
1.2.2	Trays	4	4.1.3	Operating data	39
1.2.3	Illumination unit	4	4.1.4	Analog Exit	40
			4.1.5	Collective Alarm	40
2	Delivery and Installation	5	4.2	Ordering Information	41
2.1	Delivery and Completeness Check	5	4.2.1	Configuration of the Shaking Boxes	41
2.2	Installation and Laboratory Connection	5	4.2.2	Optional Equipment	41
2.2.1	Placement at the Working Place	5	4.2.3	Accessories	42
2.2.2	Power Supply	5	4.2.4	Additional Laboratory Equipment and Consumables	47
2.3	Equipment of the Shaking Box and Optional Connections	6	4.3	Safety Instructions Laboratory Shakers and Incubation Shaking Cabinets	48
2.3.1	Shaking Drive	6	4.4	EG Declaration of Conformity	48
2.3.2	Illumination unit (optional)	6	4.5	Declaration of Decontamination	48
2.3.3	Cooling Devices	7	4.6	Information and Instructions on Disposal and Repairs	48
3	Operating Information	8	4.7	Dimension sheet CERTOMAT® BS-T	48
3.1	Safety Notes	8	4.8	Prospectus CERTOMAT® BS-T	48
3.2	Equipment, Mounting and Loading of a Tray	9			
3.2.1	Equipping the Tray	9			
3.2.2	Mounting and Dismounting the Required Tray	9			
3.3	Commissioning	10			
3.3.1	Mounting and Set-ups	10			
3.3.2	Menu Overview	11			
3.3.3	Manual Menu „MAN“	11			
3.3.4	Program menu „PROG“	15			
3.3.4.1	Program definition	16			
3.3.4.2	Example for programming	20			
3.3.4.3	Password protection for programs	21			
3.3.4.4	Cancelling of a program password	22			
3.3.4.5	Starting of a program	23			
3.3.5	Setup menu „SETUP1“	25			
3.3.6	Set-up menu „SETUP2“	31			
3.3.7	Operating of the Optional Cooling Devices	33			
3.3.7.1	Operation of an Optional External Cooling Device	33			
3.3.7.2	Operation of the Optional Internal Refrigerator	33			
3.4	Connection of External Devices	34			
3.4.1	Analog Exits (ANALOG OUT)	34			
3.4.2	RS 232 Interface	34			
3.4.3	Collective alarm	34			
3.5	Mounting of Trays and Mounting Sets	35			
3.5.1	Assembly and Loading of Finally Equipped Sets of Trays	35			
3.5.2	Mounting of Test Tube Racks	35			
3.5.3	Mounting of the Universal Mounting Set	36			
3.6	Maintenance and Trouble-Shooting	37			
3.6.1	Maintenance and Cleaning	37			
3.6.1.1	Cleaning	37			
3.6.1.2	Maintenance of the Illumination Unit	37			
3.6.2	Malfunction Correction	37			
3.6.2.1	Electrical Malfunctions	37			
3.6.2.2	Malfunctions of the Shaking Drive	37			
3.6.2.3	Malfunctions of the Optional Refrigerator	38			
3.7	Warranty Regulations and Service	38			

1 Design and Function

1.1 Design and Construction Features

The shaking box CERTOMAT® BS-T is an easy to operate, efficient and durable incubator shaking box for the universal use in biological and chemical laboratories, etc.

The basic setup of the CERTOMAT® BS-T includes the shaking drive and the integrated heater. The shaking drive consists of a motor of external rotor type with brushless drive and a noiseless, robust Poly-V belt-drive. Shaking drives with different lifting capacities are available. The lifting can be modified by our service on the job. The different shaking speeds can be adjusted continuously. The shaker can be equipped with trays of different size as well as with mounting systems of the CERTOMAT® - shaker program. For the operation within and below room temperature a cooling cycle (heat exchanger) or a refrigerating set are optionally available. For more information concerning the equipment possibilities see below, detailed information you will find in chapter 4 "Technical Data and Ordering Information". With this information you can set-up the equipment of the shaking cabinet according to your specific needs.

The box has a big cover made of Thermopane glass which allows visual control of the containers. It is easy to insert or remove the containers (Erlenmeyer flasks, containers, bottles etc.) or to change the tray. The incubation room is made of stainless steel and has a stainless steel spillage tray in the bottom. If containers break unintentionally, the spillage is sampled and cannot be released unchecked. Above this the box can be easily cleaned. The spillage tray can be easily removed of the incubation room.

1.2 Available Equipment and Accessories

1.2.1 Shaking Drives

The operating terminal and the control system are included in the front panel on right hand of the box. The shaking drive is mounted below the spillage tray bottom of the box. The drive is similar to the drive used for the bench top shaker CERTOMAT® S-II and. The drive specifications are as follows:

- Orbital motion: Ø 25 mm or 50 mm (depending on the delivered device version; can be changed to alternate amplitude in-place at customer's site by authorized service representatives of Sartorius AG.
- Shaking speeds: 40 ... 400 ¹/min
- Rotation speed accuracy: max. ±1 % of the final value

1.2.2 Trays

Trays can be delivered being finally equipped with sets of clamps for Erlenmeyer flasks of different size (type E, F), or they are available for the individual mounting as universal trays type EU or FU.

- Dimensions:

type E/EU: W x D about 420 x 420 mm¹⁾

type F/FU: W x D about 800 x 420 mm¹⁾

The universal trays can be equipped with clamps for Erlenmeyer of different size, slewable test tube racks and universal mounting sets. The universal mountings allow the setup of various vessels, such as beakers, bottles, or separation funnels, for instance.

1.2.3 Illumination unit

- The BS-T can be delivered with an optional illumination unit and includes 5 lamps of 19 W each. Each lamp can be switched off separately. The mounting of the illumination unit (see point 2.3.2) as well as the activation in the SETUP2 menu is described in point 3.3.6.

¹⁾ Trays of modified design for new locking; rebuilding of existing trays for this locking on request. Please contact the Sartorius AG customer service.

2 Delivery and Installation

2.1 Delivery and Completeness Check

The equipment and the extent of delivery of the incubation shaking box CERTOMAT® BS-T depend on the individual order. The shaking drive, and if ordered, the options „cooling cycle“ and „refrigerator“ will be preinstalled in the factory prior to delivery according to the customer's demands. Each shaking box is delivered after a thorough functional test. For the transport to the working place, the installation and prior to start-up of the unit you should check the following:

1. A suitable transport route through the building and appropriate transportation means (lifting trucks, etc.) should be available. The floor, the doors and elevators must be large enough and you should be able to move the unit without personal hazards.
2. The working place must be suited for the shaking cabinet. It must offer sufficient space for convenient placement and it must compensate the weight of the completely mounted shaking cabinet.
 - dimensions of the incubation box:
1150 x 760 x 750 mm (W x H x D)
 - weight: about 170 kg
- The BS-T can be placed on a supporting frame. This will facilitate operation.
 - Dimensions of the support frame:
1150 x 200 x 710 mm (W x H x D)
3. The working place must be able to compensate vibrations which are caused at the shaking operation. Devices which stand around may not be disturbed. Please observe the conditions for the limitation of the shaking speeds in chapter 3.

Note that the action of the shaker drive can transmit vibrations to the box and to the working place. These vibrations are caused by resonance effects which depend on the load and shaking speed. Although the drive has internal balance weights, these weights cannot fully compensate such vibrations for all operating conditions.

4. The power supply of the laboratory must be prepared according to the specifications of your CERTOMAT® BS-T, see the labels on the unit.
5. The working place should offer suitable environmental conditions, see the specifications in the supplement. If the CERTOMAT® BS-T should be applied under specific conditions or with corrosive media you will need to check the applicability of all parts and accessories. Malfunction and damages, such as corrosion by aggressive media are not subject to our guarantee.

If the working place is not applicable for the CERTOMAT® BS-T, you must not start-up the unit. Our warranty will lapse, if the shaking cabinet and its equipment are applied under inappropriate environmental conditions or with inappropriate media.

6. Please check the condition of your device and the provided accessories. All components specified in the order are premounted or are enclosed in the delivery. The delivery should be complete and all parts undamaged.
- If the CERTOMAT® BS-T has the wrong mains specifications, in case of malfunction of the unit or if any parts are damaged or are missing inform your Sartorius Stedim Biotech GmbH representative as soon as possible.

2.2 Installation and Laboratory Connection

2.2.1 Placement at the Working Place

1. Place the CERTOMAT® BS-T in the laboratory in such a way, that it is easily accessible and its operation cannot interfere with other devices placed around. The right part of the machine should be positioned in an approximate distance of 50 cm from devices or the wall. The side door is then still accessible and you must not reposition the box. A minimum wall distance of 15 cm is necessary for sufficient ventilation.
- If the shaking box contains the optional cooling cycle (heat exchanger) we recommend a distance of approximately 20 cm for the rear side.
2. The shaking box must be adjusted exactly horizontally.

Unbalanced placement of the box can cause quick wear and tear of the shaking drive, especially at heavy loads and high shaking speeds of the trays.

2.2.2 Power Supply

You will require or need to prepare a laboratory grounded socket offering 230 V, 50 Hz or 115 V, 60 Hz, depending on the country standard. Check on the type plate on the box whether it has the correct power ratings. Never connect the CERTOMAT® BS-T to power supplies of wrong ratings.

- If the plug at the power supply cable does not fit to your sockets, the service technician can attach a suitable plug during the first setup of the CERTOMAT® BS-T.

2.3 Equipment of the Shaking Box and Optional Connections

2.3.1 Shaking Drive

- The shaking drive and tray support will be mounted in the factory prior to delivery and be ready for use after assembly of the tray with the holders for the containers. If another drive system should be installed, rebuilding must only be done by authorized and qualified service personnel. Please contact your service representative or the service department of the Sartorius AG, if necessary.

2.3.2 Illumination unit (optional)

- The illumination unit (art. No. 8861463) is optionally available. We recommend to use a cooling during the operation. Without cooling the inside may become too hot due to the illumination.
- If your device has no cooling, it must be retrofitted accordingly with the models UHK (circulating air – heating – cooling device) (art. no. 8864454) or UHW (circulating air – heating – heat exchanger) (art. no. 8864560).
- Please unpack carefully the delivered illumination unit and check the delivery for damages and incompleteness.

Attention – risk of injury! You may injure at the plate edges.


Fig. 1: Mounting of the illumination unit
(Cover inside of the cabinet)

1. Built-in bar (left)
2. Fastening screws
3. Built-in bar (right)
4. Connecting cable
5. Illumination unit
6. Gas spring (right)
- 6a. Ball head with safety clip
7. Shifting square
8. Gas spring (left)
- 8a. Ball head with safety clip
9. Fastening screws

- We recommend to carry out the mounting of the illumination unit with two persons.
1. Remove the safety clips of the two ball heads (8a) of the left gas spring (8). With a rotary movement you may now remove the gas spring (8) out of the ball heads (8a).
 2. Replace the removed gas spring by the reinforced gas spring, which is part of the illumination unit.
 3. Please remove now the safety clips from the 2 ball heads (6a) of the right gas spring (6). With a rotary movement you may now remove the gas spring (6) out of the ball heads (6a).
 4. Replace the removed gas spring by the reinforced gas spring which is part of the illumination unit.
 5. Apply again the safety clips at the replaced gas springs (6/8).
 6. Completely unscrew the right built-in bar (3) in unscrewing the fastening screws (2) and release the fastening screws (9) of the left built-in bar (1) by around 1/2 rotation.
-  When inserting the illumination unit (5) please take care that the connecting cable (4) hangs down to the right side.
7. Push the illumination unit (5) from the right side laterally into the left built-in bar (1) and place it on the lower shifting square (7).
 8. Afterwards fasten the built-in bar (3), which you unscrewed before and fasten again carefully all four fastening screws (2/9).
 9. Connect the connecting cable (4) with the provided socket in the incubation room.
- After activating the option „LIGHT“ in the SETUP2-menu (see point 3.3.6) the illumination unit is ready for operation. It can be controlled manually via the front panel with key pad or programmed.

2.3.3 Cooling Devices

- If your CERTOMAT® BS-T includes the optional refrigerator (UHK), no additional connections are required. Any condensate produced during cooling operation, will be evaporated inside the cabinet and the humid air be exhausted via the ventilation of the cabinet.
 - If your CERTOMAT® BS-T includes the optional cooling cycle (heat exchanger) (UHW) the inlet and outlet connectors for the external cooling device are prepared on the rear of the cabinet. You can connect your external device via fabric hoses 15 x 2 mm. Note the labelling on the rear of the CERTOMAT® BS-T for the inlet and outlet connector.
 - A laboratory cooling water supply should deliver chloride-free water. If a FRIGOMIX® R should be applied at temperatures below +4 °C you can use commercially available cooling liquids (if questions arise on this matter, you may refer to the documentation delivered with the cooling device or contact your representative of the Sartorius Stedim Biotech GmbH).
1. Connect the supply from the cooling device to the connector on left hand (seen from the rear panel), connect the return line to the other connector. Carefully screw tight the screw nuts. Open the valves at the connectors of the CERTOMAT® BS-T by turning counterclockwise.
 2. When filling the cooling cycle for the first time, it must be carefully deaerated. The coolant should have room temperature at this time. Open the coolant supply and introduce the coolant until the cycle is completely filled and the return flow is free of bubbles.
 3. If a FRIGOMIX® R is connected, fill the coolant into the filling port of this unit, start its pump and let it run, until the return flow is free of bubbles. The filling port at the FRIGOMIX® R should also be filled. After filling the cooling cycle is ready for operation.

3 Operating Information

3.1 Safety Notes

Media used with incubation shakers can be hazardous. This is especially true for shaking cultures with living cells. If dangerous media are used only authorized personnel may have access to the unit or operate it. The operator must inform the personnel about the dangers when shaking cultures are used as well as about the necessary safety regulations.

The following safety notes must be observed:

When hazardous agents should be treated, only authorized personnel is allowed to have access to the laboratory or to the working area, respectively, and be allowed to operate the unit. It must be guaranteed that unauthorized persons do not have access to the unit.

The personnel must wear adequate working clothes and personal protection, such as gloves, protective goggles, or breath protection, for instance.

Only accessories must be used which Sartorius Stedim Biotech GmbH has released for the BS-T. Before starting operation all parts should be checked for damages, especially the glass parts. You must not use any damaged or defective (glass) equipment.

The cover with the thermopane glass is very heavy. Never let the cover fall down or bounce upward. When opening the box, carefully hold the cover and guide, until it reaches its final position. For closing, carefully guide the cover until it is closed.

Do not start the shaking drive without tray. Always switch-off the shaking drive before mounting or removing the tray or any of the containers to be incubated.

Properly fasten the tray and all mountings and holders to prevent that the containers can loosen during operation. When loading the tray you should evenly distribute the containers on the tray (i.e. balanced arrangement) and ensure that they cannot hit against each other.

Only operate the shaking box with the admissible operating values (see chapter 4 "Technical Data").

Reduce the shaking speed as far as necessary, if resonance effects occur and vibrations are transmitted to the working place or devices placed around.

At specific shaking speeds and depending on the load the action of the shaking drive can cause vibrations due to „self-resonance“ effects. Such vibrations cannot be fully compensated by the internal balance weights of the tray support and will be transmitted to the BS-T and the working place.

If media are incubated in open containers you should lower the shaking speed as far as possible to avoid unintended spillage of the media. If the media are released unintentionally, i.e. due to breaking of glass vessels, for instance, you should clean the shaking cabinet and remove the broken glass and any foreign bodies.

3.2 Equipment, Mounting and Loading of a Tray

The Sartorius Stedim Biotech GmbH offers trays for the CERTOMAT® BS-T which are finally equipped with sets of clamps for Erlenmeyer flasks of different size, or trays without assemblies. On delivery or after the first set-up of the shaking box the tray can be premounted ready for operation. If this includes all clamps and mountings required for the first use you can insert the containers which should be incubated and start operation, as shown in section 3.3. If you need to change the tray or its equipment for your application, please read the following information.

3.2.1 Equipping the Tray

- We recommend to mount all clamps and mounting sets for your containers at the tray before you place the tray in the cabinet. However, the flasks, bottles, etc. with the media (especially heavy containers) should be placed on the tray after it is attached to the tray support.

Note the weight of a tray after equipment, the tray can be very heavy. The shaking containers may disturb the handling. You can hit against the cabinet and damage the containers placed on the tray. Media can be released unintentionally and can cause personnel hazards.

1. For a universal tray without clamps, etc., place and screw tight the holders and mountings required. In order to change an existing tray equipment unscrew all parts which are not needed first.
- Clamps for Erlenmeyer flasks are available in different sizes which can be individually combined. For test tubes, bottles, beakers or separation funnels, mounting systems such as test tube racks or universal mounting systems are available. Detailed information about the available container holders and mounting systems you will find in section 3.5.
2. In order to avoid unbalances when loading the tray later on, the mountings should be fitted in that way, that the containers can be distributed evenly.

3.2.2 Mounting and Dismounting the Required Tray


Fig. 2: Mounting of the tray on the shaking table (tray shown without mounting systems and containers)

1. Place the tray onto the shaking table. Push it with the indents (1a/1b) against the rear knurled screws (1c/1d) and the clamp (2a) of the shaking table.

Please observe the weight of the tray if you have already loaded it with shaking containers.

2. Push down (2b) the front edge of the tray so that it properly fits into the front clamp (2c).
3. The borings (3a/3b) will precisely fit onto the locating pins of the tray support (3c/3d) and locate the tray.

Please take care that the locating pins (3c/3d) are placed evenly in the tray. Only then the tray is situated correctly on the shaking table.

If necessary you can adjust the locating pins as follows

1. Release the screws of the locating pins (3c/3d) with the hexagon socket screw key.
2. Eventually agitate the tray slightly. The locating pins (3c/3d) position themselves in the centre of the borings (3a/3b) of the tray.
3. Screw down the screws of the locating pins (3c/3d) carefully again after the tolerance balancing.
4. In order to improve the safe position of the tray – especially at high loadings in connection with high speed – the provided flat headed screws M5 x 8 must be screwed in the prepared borings.

3.3 Commissioning

3.3.1 Mounting and Set-ups

1. If you have changed the mounting place (e.g. cooling room <-> laboratory) please wait until the device has warmed up until room temperature before you start it. Otherwise the air humidity can condensate in the device and can lead to functional disturbances.

Following the transport of a CERTOMAT® BS-T with integrated cooling unit, the device may only be put into operation after a waiting period of an hour. This avoids danger to the cooling unit.

2. Please connect the CERTOMAT® BS-T with the mains cable to a suited socket.
3. By turning the mains switch (1) the electronics are supplied with line voltage.

After switching on the device initializes. After approx. 5 sec. the background illumination of the display as well as the green operating LED light. The manual menu "MAN" is active. The device is now ready for operation.

In pressing the ESC-key you change from the manual menu to the main menu.


Fig. 3: Frontpanel with key pad

1. mains switch
2. [+/-]-key for the set-up of process values (temperature, speed)
3. function field [F1-F5]
4. alphanumeric display
5. display field MODE
6. [ESC]-key for menu jump back
[ENTER]-key for confirmation of inputs / changes
7. Operation-LED (green)
8. [START]-key
9. ALARM-LED (red)
10. [STOP]-key

3.3.2 Menu Overview

In the menu overview you will find the menu points which are available in the main menu. In the chapters described in the following you will find the set-ups of the respectively selected menu item.


Fig. 4: Main menu (menu overview)

Operating keys

F3: SETUP-Menus: see chap. 3.3.5/3.3.6

F4: MAN manual menu: see chap. 3.3.3

F5: PROG program menu: see chap. 3.3.4

Display field:
over the MODE-field the respective "active" menu is indicated. Here: "MAIN" for main menu.

3.3.3 Manual Menu „MAN“

- In activating the key [F4] "MAN" in the main menu the „manual menu“ is activated.
- Back jump from the „manual menu“ to the main menu with [ESC]


Fig. 5: main menu


Fig. 6: MAN: manual menu

The following set-ups in the "manual menu" (fig. 5) can be carried out:

SPEED:

Set-up of speed

- After the pressing of [F1] the input field "SPEED" in speed/min. flashes.
- In pressing the [+/-] -keys the speed is adjusted.
- Confirm input with [ENTER].

TEMP:

Set-up of the temperature

- After the pressing of [F2] the input field "TEMPERATURE" flashes in °C
- In pressing the [+/-] -keys the temperature is adjusted.

Confirm input with [ENTER].

F3**TIME:****Set-up of the TIMER**

- After the pressing of [F3] the "TIMER" [TIME] can be adjusted.
- In pressing the keys [+/-] the desired operating time in [h : min] is adjusted
- Confirm input with [ENTER].
- The operating time can be set between 00h:01min and 99h:59min as well as "CONT." (continuous operation).

F4**COOLING:****Switching the cooling on and off (only at "UHK version ")**

- After pressing of [F4] the cooling is switched on / off.

Switching the cooling on is only possible at temperature set values < 38.5°C!
If a higher temperature value is entered and confirmed by pressing [ENTER] while the refrigerator is on,
a brief signal sounds and the refrigerator will be switched off.

MODE**MODE:****Display field**

- The respective "active" menu is indicated, here: "MAN" for manual menu.

- With the [F5] key you can mutually change to another display field in order to indicate all values to be set (only at option “illumination”).


Fig. 7: adjustable parameters


Fig. 8: further set-up possibilities

F2 **LIGHT:** **Switching the illumination on and off (only at option “illumination”)**

- After the pressing of [F2] the illumination is switched on / off.

MODE **MODE:** **Display field**

- The respective “active” menu is indicated, here “MAN” for manual menu.
- In pressing the key [START] you change to the operating mode “RUN”.
- The previously adjusted values for:

Speed	:	SPEED
Temperature	:	TEMP
Lighting	:	LIGHT
Cooling	:	COOLING
Timer	:	TIME

are activated now.

- With the key [F5] you can alternately change to another display field in order to indicate all actual values.


Fig. 9: Operating mode RUN

< F5 >


Fig. 10: Operating mode RUN

- The “TIMER“ runs. The respective rest operating time is indicated.

In the operating mode the ACTUAL values are indicated.
Set values cannot be changed here!

- In pressing the function keys [SPEED], [TEMP], [COOLING] and [LIGHT] a little bit longer the parameter flashes and the SET-value is indicated. After the release of the function keys the ACTUAL-values are indicated again.
- The ACTUAL value display of the temperature refers to the air temperature in the incubation room.
- In pressing the key [STOP] you jump back to the manual menu “MAN“.
- Heating, cooling, light and shaker are switched off.
- A new set-up can be effected.
- After the adjusted operating time is finished, you automatically jump back to the manual menu “MAN“.

Once key [STOP] has been pressed, the refrigerator will continue to operate for up to 3 minutes. The refrigerator can only be restarted after 3 minutes.

3.3.4 Program menu "PROG"

- In activating the key [F5] "PROG" the „program menu“ is activated. "PROG1" flashes!
- Jump back to the main menu with [ESC]


Fig. 11: main menu


Fig. 12: PROG1-4: program menu

- With [F1][F5] the respective program is selected. After the selection the selected program flashes in the display.
- In the „MODE“-field the setup program appears automatically, here "PROG1".
- The programming according to an example is explained in point 3.3.4.2 "Example for programming".

3.3.4.1 Program definition

- After the pressing of the [ENTER] key (alternately press the flashing function key for approx. 2 sec.) you change to the program definitions "PRSTEP . .STEP1 . .STEP4".
- In this menu the different "STEPS" can be selected within the program "PROG1.....PROG4".
- At programs protected by passwords the respective password must be input for the editing/modification of present programs.


Fig. 13: PROG1-4: program menu


Fig. 14: program definition

- In fig. 14 (program definition) is indicated, whether the PRSTEP is "ON" or "OFF". Furthermore the step duration is indicated.

Selection of the different "STEPS" (fig. 14):

F1 for: **PRSTEP**

- After the pressing of [F1] the set-ups for "PRSTEP" can be effected.

F2 for: **STEP1**

- After the pressing of [F2] the set-ups for "STEP1" can be effected.

F3 for: **STEP2**

- After the pressing of [F3] the set-ups for "STEP2" can be effected.

F4 for: **STEP3**

- After the pressing of [F4] the set-ups for "STEP3" can be effected.

F5 for: **STEP4**

- After the pressing of [F5] the set-ups for "STEP4" can be effected.

MODE **MODE:** **Display field**

The respective "active" PROG is indicated, here "PROG1".

- With the keys [F1 to F5] the adjustable parameters of the steps can be selected.


Fig. 15: program definition



Fig. 16: adjustable parameters

For each “STEP“ the following parameters can be adjusted (fig. 16):

F1

SPEED:

Set-up of the speed

- After the pressing of [F1] the input field “SPEED“ in speed/min. flashes.
- In pressing the keys [+/-] the speed is adjusted.
- Confirm input with [ENTER].

F2

TEMP:

Set-up of the temperature

- After the pressing of [F2] the input field “TEMPERATURE“ in °C flashes.
- In pressing the keys [+/-] the temperature is adjusted.
- Confirm input with [ENTER].

F3

TIME:

Set-up of the STEP-duration [h : min]

- After the pressing of [F3] the input field “TIME“ flashes.
- In pressing the keys [+/-] the time interval (duration) is adjusted.
- Confirm input with [ENTER].
- At the STEP duration values between 00h:01min (immediate change to the next step) and 99h:59min as well as “CONT.“ (continuous operation until the manual break with STOP) can be adjusted.

You cannot determine a firm operating time for the PRSTEP. It can only be switched “ON“ or “OFF“!

- PRSTEP in operating mode: the conditions selected in the PRSTEP from the pressing of the [START] key up to the actual program start [STEP1] are activated.
- PRSTEP not operating: the shaker, tempering, cooling and light are deactivated until the program start [STEP1].

F4**COOLING:****Switching the cooling on and off (only at "UHK version")**

- After the pressing of [F4] the cooling is switched on / off.

Switching the cooling on is only possible at temperature set values <38.5 °C!
 If a higher temperature value is entered and confirmed by pressing [ENTER] while the refrigerator is on,
 a brief signal sounds and the refrigerator will be switched off.

MODE**MODE:****Display field**

- The respective "active" STEP is indicated.
- With the [F5] key you can alternately change to another display field in order to set up further parameters (only at option „illumination“).


Fig. 17: adjustable parameters


< **F5** >

Fig. 18: further set-up possibilities

F2**LIGHT:****Switching the illumination on and off
(only at "option illumination")**

- After the pressing of [F2] the illumination is switched on / off.

MODE**MODE:****Display field**

- The respective "active" STEP is indicated, here "STEP1"

- In pressing the key [ESC] you jump back to the previous menus (until the main menu)


Fig. 19: set-up possibilities within a „STEP“


Fig. 20: program definition


Fig. 21: program menu


Fig. 22: main menu

3.3.4.2 Example for programming

- Altogether five different programs can be memorized in the device. Each program consists of one "PRSTEP" and four further steps which can be processed one after the other. For each step the values for cooling, light, speed and temperature can be set independently from one another. Each step has a defined duration time.
- If only one or two steps are required the remaining steps can be defined with a duration time 00:00 so that the program is directly stopped afterwards. A repeat function for endless loops is available.
- The "PRSTEP" can only be switched on or off (no operating time). It defines the conditions (cooling, light, speed and temperature) from the pressing of the key [START] until the program start (STEP1).

- At "PRSTEP" switched off, the cooling, light, speed and temperature are switched off until the program start.
- The main purpose of use for the "PRSTEP" is the date/time exactness for the provision of cultures, e. g. after the weekend or the next morning.

Advantages:

- You have to define the program only once with operating times and you can start it afterwards at any time. The time of the program end remains always the same.
- Comfortable start of night or weekend programs during the day or before the weekend.
- Conditions are reproducible before the start of the program because these are held by the device according to the settings in the "PRSTEP".
- Automatic program start without the personnel being present.

Example:

On Monday morning at 9:00 o'clock any cultures have to be ready. The program (STEP1...STEP4) for the production of these cultures needs altogether 38 h. So the program must be started manually without use of the "PRSTEP" on Saturday evening at 7:00 o'clock.

The "PRSTEP" makes it possible to choose any time for the start and to keep the material cooled under the desired conditions.

PROG1: PRSTEP: ON at 20 °C, speed 0, light off, cooling on

Step 1: 7 h at 30 °C, speed 200, light on, cooling off

Step 2: 10 h at 40 °C, speed 100, light off, cooling off

Step 3: 20 h at 10 °C, speed 60, light on, cooling on

Step 4: 1 h at 10 °C, speed 200, light on, cooling on

Function	PRSTEP	STEP1	STEP2	STEP3	STEP4
Cooling	on	off	off	on	on
Light	off	on	off	on	on
Speed [rpm]	0	200	100	60	200
Temperature [°C]	20	30	40	10	10
TIME (operating time)	ON	7 h	10 h	20 h	1h
Day	Saturday	Sunday	Sunday	Monday	Monday
Time	19:00	02:00	12:00	08:00	09:00
	START PRSTEP	Program START			Program end

Fig. 23: Example for programming

3.3.4.3 Password protection for programs

- The “password protection” (code) for programs serves for the protection of programs which have been once established for this purpose against unauthorized access. The content of the program cannot be modified without knowledge of the password. If requested also the “starting” and “stopping” of the programs can be protected by this password.
- After the pressing of the [ENTER] key you change to the menu where you input the “password” (code) for the respective active program.
- For safety reasons the “password” must be input twice if the password is determined.


Fig. 24: program definition of “PROG1”

ENTER >


Fig. 25: determination of password for “PROG1”

- In pressing the function keys [F1] ... [F5], [+], [-], [START] or [STOP] as admissible signs the “password” (code) (visible “★” on the display) is given in.
- The “password” (code) must have at least 4 and at most 7 places.
- After the input of the password the “password” (code) must be confirmed with the [ENTER] key. Then you change automatically to the confirmation of the password which has been input.
- After the repeated input of the password the “password” (Code) must be confirmed once again with the [ENTER] key in order to secure the “password input”. At the same time the back jump to the “program definition menu” is effected.
- In pressing the [ESC] key the “password determination” can be stopped. All “password inputs” which have been effected until this time become invalid.

3.3.4.4 Cancelling of a program password

- Select in the program definition menu the program in which you want to cancel the password (password knowledge is necessary).
- In pressing the [ENTER] key you change from the program definition menu to the menu for the password determination.


Fig. 26 program definition of "PROG1"

ENTER >


Fig. 27: determination of password for "PROG1"

- Deleting / deactivating of the "password" (code) is carried out with the [ENTER] key as first password sign. You change then to the confirmation of the deletion / deactivation.
- The repeated input of [ENTER] causes the deletion of the existing "password" (code) and the automatic back jump to the program definition menu is effected.
- In pressing the [ESC] key the deletion / deactivation of the "password" (code) can be stopped.

3.3.4.5 Starting of a program

- In pressing the [START] key afterwards you change to the setup mode of the program selected before, here program 1 "PROG1".


Fig. 28: PROG1-4: program menu


Fig. 29: start menu

The following parameters in the start menu (fig. 29) can be adjusted:

F1

TIME:

Set-up of the time (start time of the program)

- After the pressing of [F1] the input field "TIME" flashes.
- In pressing the [+/-]-key the time is adjusted.
- Confirm input with [ENTER].

F2

DAY:

Set-up of the date (start date of the program)


- After the pressing of [F2] the input field "DAY" flashes.
- In pressing the [+/-] keys the date is adjusted.
- The current date and the current time are respectively preadjusted.
- Confirm input with [ENTER].

F3

REPEAT:

Set-up of repeat

- After the pressing of [F3] the automatic program repeat is activated / deactivated.
- At activated program repeat the "STEP 1" to "STEP 5" of the program is processed in an endless loop.
- The "PRSTEP" is cycled only at the first time.

- The respective “active” program is indicated, here “PROG1”.
- In pressing the [START] key afterwards you change to the operating mode “RUN”. The program activated before is processed.
- With the [F5] key you can alternately change to another display field in order to indicate all actual values!


Fig. 30: RUN: operating mode





Fig. 31: RUN: operating mode

- At the TIME-display the rest operating time of a “STEP” is indicated. During the “PRSTEP” processing no time is indicated (display “PRSTEP”).
- In the operating mode the ACTUAL VALUES are indicated and the started program runs automatically.
- In pressing the respective function keys the SET VALUE is indicated.
- After having released the function key the ACTUAL VALUES are indicated again.
- In pressing the [STOP] key you jump back to the program menu “PROG”. Heating, cooling, light and shaking drive are switched off. The program processing is stopped.

It is indicated during the program processing which program and which “STEP” are being processed at the moment, e. g.:

- For the programs “PROG2-PROG5” the set-ups are used analogously to the description for program 1 “PROG1”.

3.3.5 Setup menu „SETUP1“

- In pressing the [F3] key you change to the “setup1 menu”.
- All set-ups and modifications are memorized internally and remain active also after voltage loss or switching off and on!


Fig. 32: main menu


Fig. 33: setup1 menu

Display contrast

- Can be modified in pressing the keys [+] and [-].

In pressing the key [-] the contrast is so much reduced that no information is visible at the display!

The following parameters in the “setup1 menu“ (fig. 33) can be adjusted:

F1

LANG.:

Set-up of the menu language

- After the pressing of the [F1] key you come to the menu in which you adjust the menu language (German, English, French).
- An arrow marks the respectively active language.
- In pressing the keys [+] / [-] the arrow can be moved.
- In pressing the [ENTER] key the respectively adjusted language is confirmed.
- In pressing the [ESC] key you change again to the “setup1-menu“. The language set-up before is now shown in the menus.

F2**DIAG:****Diagnosis of the alarm messages which occurred**

The alarms which occurred are indicated on a list. The following information is separately made for disposal for each error type:

Error type	e. g. speed or temperature alarm
Error counter	Number of the alarms which occurred since the last acknowledgment
Date / time / alarm start	Date and time of the alarm start (only of the alarms which occurred last)
Date / time / alarm end	Date and time of the alarm end (only of the alarms which occurred last)

- An arrow marks the respective active list entry.
- In pressing the keys [+/-] the arrow can be moved.
- In pressing the [ENTER] key the respectively active list entry is cancelled. Afterwards the alarm information for this kind of error is not available any more.
- If no alarm occurred or if all list entries were cancelled [„.....<ESC>....“] appears on the display. The diagnosis menu can only be quit with [ESC].

F3**ALERTS:****Alarm definition**

In a list the alarm parameters for temperature, speed, door contact and voltage loss can be adjusted.

- An arrow marks the respectively active parameter.
- In pressing the keys [+/-] the arrow can be displaced.
- In pressing the [ENTER] key the YES / NO definitions are directly changed or the numeric values can be released for modification (value flashes).
- In pressing the keys [+/-] the respectively flashing value can be set.
- In pressing the [ENTER] key the adjusted value is confirmed.
- In pressing the [ESC]-key you jump back to the SETUP1-menu.

The following alarm parameters can be set up:

Temperature 1	ON	The temperature alarm is active
	OFF	The temperature alarm is inactive
Signal	ON	Acoustic signal active in case of an alarm
	OFF	Acoustic signal inactive in case of an alarm
Time OK	40 [s]	Duration of time, during which the set temperature must have been at least reached after a temperature error, before the next error can be registered (here 40 seconds).
Time ERR	40 [s]	Duration of time, during which the temperature error must at least occur, before the error is indicated (here 40 seconds).
min.	3.0 [°C]	Permitted temperature deviation downwards, where no error is indicated (here 3,0°C)
max.	3.0 [°C]	Permitted temperature deviation upwards, where no error is indicated (here 3,0°C)
Door contact		
	ON	The door contact alarm is active
	OFF	The door contact alarm is inactive
Signal	ON	Acoustic signal is active if the alarm occurs
	OFF	Acoustic signal is inactive if the alarm occurs
Time OK	5 [s]	Duration of time, during which the plexiglass cover must be at least closed, before the error is registered again (here 5 seconds)
Time ERR	2 [s]	Duration of time, during which the plexiglass cover must be at least opened, before the error is indicated (here 5 seconds)
Speed		
	ON	The speed alarm is active
	OFF	The speed alarm is inactive
Signal	ON	Acoustic signal active at alarm message
	OFF	Acoustic signal inactive at alarm message
Time OK	20 [s]	Duration of time, during which the set speed after a speed error must have been at least reached again, before the next error can be registered (here 20 seconds)
Time ERR	20 [s]	Duration of time, during which the speed error must at least occur, before the error is indicated (here 20 seconds)
min.	10 [U/min]	Permitted speed deviation downwards, where no error is indicated (here 10 U/min.)
max.	10 [U/min]	Permitted speed deviation upwards, where no error is indicated (here 10 U/min.)
Voltage loss		
	ON	The voltage loss alarm is active
	OFF	The voltage loss alarm is inactive
Signal	ON	Acoustic signal is active if the alarm occurs
	OFF	Acoustic signal is inactive if the alarm occurs

- In the case of an alarm there is an optical display via alarm-LED or an acoustic message via buzzer (if active). At the same time the alarm is taken over in the diagnosis and can be read-out after the program end in the point SETUP / DIAG.
- A message to outside can be sent over a contact without potential via the SUB-D socket (see point 3.4.3).
- In pressing the [ENTER] key the optical and acoustic messages are acknowledged as well as the contact without potential is set back.
- In opening the doors the shaking drive is stopped and the heating is switched off, independent of the fact, if the door contact alarm is ON or OFF. If the cover is closed again, the shaking drive and the heating are operating again automatically.
- In the case of voltage loss during the program processing or in the operation mode "RUN" the time which has already been processed is saved. If there is voltage again, the device starts automatically. The remaining time of the program is processed. This happens independent of the fact if the voltage loss alarm is active or not.
- The alarm definition for the voltage loss serves only for the activation of the optical/acoustic alarm message or the entry in the diagnosis after the voltage has returned.

F4 TIME: Set-up of time and date

- In pressing the key [F4] you can change to the display field “set-up of time / date”.


Fig. 34: setup menu



Fig. 35: set- up time/date

The following parameters for the set-up of “time and date” (fig. 35) can be adjusted:

F1 DAY: Set-up of the week day

- After pressing of [F1] the input field “DAY“ flashes.
- In pressing the keys [+/-] the desired week day can be set up.
- Confirm input with [ENTER].

F2 MONTH: Set-up of the month

- After the pressing of [F2] the input field “MONTH“ flashes.
- In pressing the keys [+/-] the desired month can be set up.
- Confirm input with [ENTER].

F3 YEAR: Set-up of the year

- After pressing of [F3] the input field “YEAR“ flashes.
- In pressing the keys [+/-] the desired year can be set up.
- Confirm input with [ENTER].

F4

HOUR:

Set-up of hours

- After the pressing of [F4] the input field "HOUR" flashes.
- In pressing the keys [+/-] the desired hour can be set up.
- Confirm input with [ENTER].

F5

MINUTE:

Set- up of the minutes

- After the pressing of [F5] the input field "MINUTE" flashes.
- In pressing the keys [+/-] the desired minute can be set up.
- Confirm input with [ENTER].

MODE

MODE:

Display field

- The respective "active" menu is indicated, here "SETUP1".

3.3.6 Set-up menu “SETUP2”

- In pressing the [F5] key you change to the “SETUP2-menu” for the input of the device parameters.
- All set-ups and modifications are memorized internally and remain active also after voltage loss or switching off and on!


Fig. 36: setup1 menu


Fig. 37: setup2 menu

The following “device paramters“ (fig. 37) can be adjusted:

F1

DEVICE:

Activation of different device options

Device options can be activated or deactivated in a list as well as device information can be called.

- An arrow marks the respectively activated parameter.
- In pressing the keys [+/-] the arrow can be moved.
- In pressing the [ENTER] key the YES / NO-definitions are directly modified or numeric values are released for modification (value flashes).
- In pressing the keys [+/-] the respectively flashing value can be adjusted.
- In pressing the key [ENTER] the adjusted value is confirmed.
- In pressing the [ESC] key you jump back to the SETUP2 menu.

Temp.Offset	0,0	Actual value correction for the 1.PT100 (for temperature regulation). Calibrating values can be adjusted. The digit sign can be modified by pressing the [F1 ... F5] keys.
Light	Yes	Illumination active
	No	Illumination inactive
August 2002	111	Software version
Ser No.	39	Hardware serial number (here no. 39)

- In pressing the [F2 und F3] keys you change to menus where you input internal parameters.


Fig. 38: setup2 menu


Fig. 39: password input menu

These menus are reserved for service purposes and are protected by passwords.

Modifications of these parameters can only be carried out by authorized personnel.

Improper modifications may destroy the unit.

- In pressing the [F4] key a little bit longer the device is reset.

All set-ups, programs, passwords etc. will definitely get lost!

3.3.7 Operating of the Optional Cooling Devices

3.3.7.1 Operation of an Optional External Cooling Device

- If the shaking cabinet is equipped with the cooling cycle (heat exchanger) you must activate the external cooling device or the cooling water supply manually or switch it off. Please make sure that the heating of the BS-T and the external heating work against each other.
 - The cooling water temperature and flow rate determine the temperature in the incubation room of the CERTOMAT® BS-T. The possible incubation temperature is within 10 ... 15 °C above the cooling water temperature. If a cooling device is connected, you can select the temperature within the specifications of this device, i.e. -10 ... +40 °C for a FRIGOMIX® R.
 - For a temperature of below +10 °C we recommend to use a commercially available coolant. Otherwise ice may form in the device, which may affect the cooling function.
1. Check the connections of the cooling water supply or the connections of the external cooling device (supply and return flow) and check whether the circuit is completely filled.
 2. Open the manual valves on the rear connectors of the CERTOMAT® BS-T and switch-on the external cooling device (or coolant supply), if cooling is required. Ensure that the coolant can circulate properly. The tubings must not be bent or clamped off unintentionally.
- If bubbles occur in the return flow, the cooling cycle should be refilled.
3. If a FRIGOMIX® R is applied, adjust the set temperature for cooling. Turn the right rotary knob while pressing the red push button on the left.
 4. If the set temperature of the CERTOMAT® BS-T is > 5 °C above the room temperature or higher, you should switch off the external cooling device and/or close the supplies.
- In addition you should close the valves at the connectors on the rear of the CERTOMAT® BS-T. For this turn the valves clockwise. Then you can disconnect the supply tubings without emptying the heat exchanger inside the CERTOMAT® BS-T.

3.3.7.2 Operation of the Optional Internal Refrigerator

- If the CERTOMAT® BS-T is equipped with the optional refrigerator, you need not care for condensate produced during cooling operation. Existing condensate will be evaporated internally in the cabinet.
- To avoid damages of the refrigerator an integrated overheat protection permanently switches off the unit if the temperature in the incubation room exceeds +38.5 °C.

3.4 Connection of External Devices

At the rear side of the device the CERTOMAT® BS-T has two Sub-D sockets with 9 poles for analog signal exits, RS 232 interface, as well as a collective alarm contact.

3.4.1 Analog Exits (ANALOG OUT)

- For the external registration of speeds and temperatures, e. g. with a recorder, the CERTOMAT® BS-T has analog signal exits (0 ... 10 V).

The conditions are as follows:

- Speed: 0 V corresponds with 0 min⁻¹
10 V corresponds with 400 min⁻¹
- Temperature: 0 V corresponds with 0°C
10 V corresponds with 100°C

The moving load at the voltage entries must not be lower than 10 kΩ.

- Optionally current signals (0 ... 20 mA) or (4 ... 20 mA) can be made for disposal.

The moving load at the current entries must not be higher than 500 kΩ.

3.4.2 RS 232 Interface

- The Sub-D socket with 9 poles "RS 232" is currently a Sartorius AG service interface.

3.4.3 Collective alarm

- A contact without potential is made available on PIN 4/9 via a Sub-D-Socket "ANALOG OUT". Alarm situations can be reported externally by help of this contact.

Contact loading capacity: 230 V AC / 0.5 A (Ohm load).


Fig. 40: Device rear side

- RS232 interface
- Socket "Analog Out" with 9 poles for recorder connection etc.
 - Pin 1 = signal + / temperature
 - Pin 2 = Gnd / temperature
 - Pin 6 = signal + / speed
 - Pin 7 = Gnd / speed
 - Pin 4 = collective alarm
 - Pin 9 = collective alarm
- Supply entry socket with supply main and shock-proof plug

3.5 Mounting of Trays and Mounting Sets

3.5.1 Assembly and Loading of Finally Equipped Sets of Trays


Fig. 41: Tray type F with clamps for Erlenmeyer flasks

The trays type E and F are available as completely equipped tray. Detailed information you will find in the "Ordering Information" in the supplement.

1. Place the tray in the BS-T and fasten it as explained in point 3.2.2 (Mounting and Dismounting a Tray).
2. Insert the Erlenmeyer flasks in the clamps (you can also use bottles, beakers etc., which fit into the clamps).

Please take care of the distribution of the loads and distribute them evenly starting from the centre of the tray outward.

Before starting the shaker check the proper seat of all containers. They must not get loose and flung off during operation of the shaker.

3.5.2 Mounting of Test Tube Racks


Fig. 42: Test tube rack assembly

1. Loosen the clamping screw of the test tube rack and turn the test tube holder, until the screw-holes for attachment of the rack to the tray are accessible. If necessary, you may completely remove the holder from the rack.
2. Place the test tube rack onto the universal tray and tighten with the mounting screws. Consider the number of racks required for your application and distribute the racks evenly on the tray starting from the centre of the tray outward.
3. As far as being disassembled, mount again the test tube holders to their racks. Turn the holders into their working position (usually vertical). Locate with their clamping screw.

3.5.3 Mounting of the Universal Mounting Set

The universal mounting set is available for the trays of size EU and FU. Enclosed you will find further information about the component parts. Basic parts are two side frames which are adapted to the corresponding tray and the rods for connection of the side frames:

- Side part set type B-2, art.-no. 8854238, for tray EU
- Side part set type B-3, art.-no. 8854243, for tray FU

Available universal clamping rods and mounting sets for separation funnels are as follows

- Universal clamping rods type U-2, art. no. 8854289, for side part set B-2/B-3
- Mounting set type S-1 for separating funnels 50 and 100 ml, art. no. 8854262, with clamps and retaining springs and 1 universal clamping rod type U-2.
- Mounting set type S-1 for separating funnels 250, 500 and 1000 ml, cat.-no. 8854270, with clamps and retaining springs and 1 universal clamping rod type U-2

1. Place the side parts (1) onto the tray and screw them onto the tray with the mounting screws as shown at pos. (1a). For each side part 4 screws are supplied to provide proper tightening. Then mount connecting rods as shown at pos. (1b).
2. For a mounting system with universal clamping rods attach the clamping rods at the side frame rods as shown at pos. 3. You'll need at least two rods and you can add as much as required for placement of your vessels.
3. Place the shaking containers between the clamping rods. Distribute the containers evenly on the tray, starting from the centre outward. You can loosen the clamping screws 3a/3b of the rods to readjust the rod and properly locate the flasks, vessels and beakers, etc.. Then carefully tighten the clamping screws 3a/b so that the containers do not loosen during operation.
4. For the mounting set for separating funnels attach one universal clamping rod as shown above. It will serve as support. Then mount the retaining springs (2a) to the clamping rod for separation funnels (2) and attach this part to the side part assembly. Depending on the size of the separation funnels you can mount up to 5 retaining springs to the clamping rod (2a).
5. Place the separating funnels onto the mounting set and locate properly. You can readjust the clamping rods as required using their clamping screws attached.


Fig. 43: Assembly of universal rods and mounting sets for separating funnels

3.6 Maintenance and Trouble-Shooting

3.6.1 Maintenance and Cleaning

Except for necessary cleaning, the replacement of defective fuses and checking of the shaking drive belt the CERTOMAT® BS-T is maintenance-free. If the optional illumination unit is installed, you can replace defective lamps. All bearings of the shaking drive are lubricated for their lifetime. The drive motor is a service-free brushless motor of external rotor type.

Disconnect the mains cable from laboratory grounded socket before cleaning and/or maintenance and service. Unintended starting of the shaking drive during cleaning and maintenance can cause personnel hazards.

- Maintenance and service of the shaking drive, such as the replacement of the drive belt, service of the electrical and electronic equipment or of the optional cooling devices must only be done by qualified and authorized service personnel.

3.6.1.1 Cleaning

1. The interior incubation room consists of stainless steel. For cleaning, we recommend usual household cleaning agents or alcohol. Do not use aggressive cleaning agents. Chlorous (hypochlorite) agents, for instance, can cause corrosion.
2. After incubation of living microorganisms and cells the CERTOMAT® BS-T may need to be disinfected. For disinfection we recommend disinfectants based on alcohols, such as Sartorius Stedim Biotech GmbH Meliseptol®, see ordering information in the supplement. After cleaning and disinfection with such agents you should thoroughly dry and ventilate the incubation room.

Remove broken glass, spilled liquids and impurities immediately after they occur.

3. The spillage tray integrated in the bottom of the shaking box can be removed for cleaning purposes.
4. Take care that no foreign bodies can enter the air-circulation system. Remove dust from the outer covering grids of the built-in ventilator and of the refrigerator regularly.

3.6.1.2 Maintenance of the Illumination Unit

- Maintenance is confined to the checking for and replacement of defective fluorescent lamps and device fuses. If other malfunctions occur, please contact your service representative of Sartorius AG.
1. Switch-off the illumination unit and disconnect from its mains supply.

2. Pull the defective fluorescent tube out of its holder. Insert the new tube. On delivery commercially available fluorescent tubes „L18W/77 (Fluora)“, or fluorescent tubes according to the customer's specifications are installed.
3. Replace defective fuses, if necessary. The fuses are placed in a fuse assembly behind the door at the right hand side panel. Note to use the same type/size for replacement, as delivered from factory side.

3.6.2 Malfunction Correction

3.6.2.1 Electrical Malfunctions

1. If the BS-T does not start, check the mains supply. The mains cable should be properly connected and the laboratory mains are operative. If other devices in the laboratory are malfunctioning at this time, a mains failure (such as voltage fluctuations) can be the cause.
 2. Check the fuses. They are located in an insert inside the cabinet, behind the right side panel. A burned-off body or a melted wire will indicate a defective fuse. For replacement note the correct voltage/current of the fuses.
- If the BS-T cannot be started although the mains supply is operative and the fuses are in good condition, contact your service representative of the Sartorius AG. You will find the telephone number on the next page.

3.6.2.2 Malfunctions of the Shaking Drive

- If the shaking drive runs noisy or irregularly the belt can be dirty. If the tray support does not turn although the shaking drive is working (you can hear the motor run), the drive belt may be damaged. If the motor is inoperative although the power supply and all necessary adjustments are available, the motor may be damaged.
1. For checking the belt drive the tray and tray support must be disassembled. Unscrew the 4 mounting screws in the centre of the tray support and remove the table. You can check the belt through the opening in the cover of the drive.
 2. If you mount the shaking table again, please carefully tighten the screws again.
- If the drive belt is damaged (torn) or the motor is suspected to be defective, please contact your service representative of the Sartorius AG, see below.

3.6.2.3 Malfunctions of the Optional Refrigerator

- The refrigerator includes an internal winding protection for overheat detection. This device cuts off the power supply at overheating of the winding. The refrigerator will automatically be enabled for restart after cooling down of the winding. This may last up to 1 hour.
- If the overheat protections cuts-off the refrigerator often or repeatedly, you should contact your service representative of the Sartorius AG.
- Make sure that the ventilation grids did not sink in and the device has a sufficient wall distance.

3.7 Warranty Regulations and Service

The Sartorius Stedim Biotech GmbH warrants its products according to the terms specified in the „General Terms and Conditions of Business“ and unless other terms were agreed upon in writing. Date of reference is the date of delivery.

- The guarantee includes defects and malfunctions due to construction, production and material faults, however it does **not** include defects due to wrong operation and improper treatment and use. Furthermore it does not apply for parts which became damaged due to normal wear and for consumer materials.
- The guarantee will lapse if the user or unauthorized third persons modify the device technically or if they use equipment, components and accessories, which have not been passed by Sartorius Stedim Biotech GmbH for the application in connection with the shaking cabinet.

Repairs can be carried out by the Sartorius AG service representative or by authorized personnel of a customer-owned workshop. Please inform your local Sartorius AG sales representative or contact:

Sartorius AG
Servicezentrum Nord
Weender Landstraße 94-108
37075 Göttingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Telefax +49.551.308.3730

- Only spare parts passed by Sartorius Stedim Biotech GmbH for the CERTOMAT® BS-T must be used.
- Defective units can be returned to the Sartorius AG.

Those apparatus sent to us for repair have to be in hygienic perfect condition, clean and well packed. If parts were polluted by media and media components, they have to be cleaned, decontaminated, disinfected or sterilised, according to the valid safety rules of this particular application, e.g. for chemical or biological safety.

The sender has to prove that he observed all safety guidelines. For this purpose you have to fill in the enclosed decontamination declaration and attach it to the apparatus. No apparatus will be repaired without declaration of decontamination or description of the measures taken.

Transport damages or subsequent cleaning or disinfection, if required, will be borne by the sender.

4 Technical Data, Ordering Information

4.1 Technical Data

4.1.1 Mechanical Design

Shaking box	outside housing 1150 x 760 x 750 mm (W x H x D)
Weight (with equipment)	about 170 kg
Tray support	440 x 420 mm (W x D)
Materials	outside cabinet: sheet steel case with acidproof varnish; incubation room inside: stainless steel case
Trays, type/size	Type E/EU (420 x 420 mm) ¹⁾ ; type F/FU (800 x 420 mm) ¹⁾
Load	depending on tray and mounting system of the Sartorius Stedim Biotech GmbH accessory program

4.1.2 Mains Supply

Laboratory connector	laboratory grounded socket
Mains specifications	230 V (±5 %) / 50 Hz or 115 V (±5 %) / 60 Hz
Power consumption	about 1.2 kW including cooling
Illumination	about 90 W (5 x 18 W)
Fuses	T6.3A for 230 V, T10A for 115 V
Radio interfer. suppress	according to EN 55011 and EN 55014

4.1.3 Operating data

Orbital motion	Ø 25 mm or 50 mm, according to shaking drive	
Shaking speed	40 ... 400 U/min	
Rotation speed accuracy	max. ±1% (from max. value))	
Working temperature	+5 °C over room temperature ... +70 °C	
Incubation temperature	Heating version	RT +8 °C ... +70 °C
	Heating and cooling version	RT -10 °C ... +70 °C
Ambient temperature	+10 ... +60 °C (+10 ... + 30 °C at UHK version)	
Ambient humidity	10 ... 60 %; incubation room 10 ... 80 %, not condensing	

¹⁾ Trays of modified design for new locking; rebuilding of existing trays for this locking on request. Please contact the Sartorius AG customer service.

4.1.4 Analog Exit

Socket, 9-poles Sub-D 0 ... 10 V (min. moving load = 10 k Ω)
(Pin configuration see point 3.4.1, fig. 40)
0 ... 20 mA or. 4 ... 20 mA (optionally)
(max. moving load = 500 Ω)

Exactness ± 0.3 V

4.1.5 Collective Alarm

Collective alarm Contact without potential (closer), max. 230 VAC (0.5 A Ohm load) via SUB-D socket
"analog out" pin 4/9

(We reserve the right to carry out dimension and construction modifications)

4.2 Ordering Information

4.2.1 Configuration of the Shaking Boxes

Cat.-no.	Ordering data, features and specifications
	Inkubation shaker CERTOMAT® BS-T / 25 mm Version 230 V / 50 Hz

8865426	Heating version (UH)
8865620	Heating and cooling version (UHK)

Inkubation shaker CERTOMAT® BS-T / 50 mm
Version 230 V / 50 Hz

8865523	Heating version (UH)
8865728	Heating and cooling version (UHK)

4.2.2 Optional Equipment

Cat.-no.	Ordering data, features and specifications
----------	--

8864489	Support frame Squared frame construction with sheet metal panel. shock-resistant stove enamel finish. Dimensions W x H x D: 1150 x 200 x 710 mm. 4 height adjustable stands
---------	--

8861463	Illumination unit for CERTOMAT® BS-1 / BS-T consisting of five fluorescent lamps à 18 W can be switched of separately, programmable Version 230 V / 50 Hz – in connection with cooling UHK
On request	Version 115 V / 60 Hz – can be integrated afterwards


8864560	Kühlregister für BS-1 / BS-T Heat exchanger (UHW) for the operation with external cooling, e. g. with FRIGOMIX® R or external cooling water supply to be delivered ex factory
---------	---

4.2.3 Accessories

Cat.-no.	Features and Specifications
Trays with fixed clamps (with clamps for Erlenmeyer flasks)	
Tray E 420 x 420 mm – clamps made of stainless steel	
8853533	Tray incl. 39 clamps for Erlenmeyer flasks 100 ml
8853568	Tray incl. 20 clamps for Erlenmeyer flasks 250 ml
8853584	Tray incl. 14 clamps for Erlenmeyer flasks 500 ml
8853606	Tray incl. 9 clamps for Erlenmeyer flasks 1.000 ml
Tray E 420 x 420 mm – plastic clamps	
8853688	Tray incl. 39 plastic clamps for Erlenmeyer flasks 100 ml.
8853666	Tray incl. 19 plastic clamps for Erlenmeyer flasks 250 ml.
8853677	Tray incl. 14 plastic clamps for Erlenmeyer flasks 500 ml.
Tray F 800 x 420 mm – clamps made of stainless steel	
8853738	Tray incl. 74 clamps for Erlenmeyer flasks 100 ml
8853762	Tray incl. 40 clamps for Erlenmeyer flasks 250 ml
8853789	Tray incl. 26 clamps for Erlenmeyer flasks 500 ml
885 3800	Tray incl. 15 clamps for Erlenmeyer flasks 1.000 ml
Universal - Trays	
8853002	type EU (420 x 420 mm)
8853037	type FU (800 x 420 mm)
Sticky tape for universal trays	
8860416	Sticky tape "Premium" for universal trays, 30 x 1 mm, roll of 10 m, durable quality
8864497	Sticky tape for universal trays, size 30 x 1 mm, roll of 50 m for fixing of Erlenmeyer flasks, beakers, etc., only for use at low shaking speeds
Anti-skid layer for universal trays	
8864470	Anti-skid layer size 380 x 420 mm for placement of Erlenmeyer flasks, beakers, etc., only for use at low shaking speeds and temperatures up to 50 °C

Test Tube Racks

Fig. 44:
Test tube rack

8853134	test tube rack for 64 test tubes Ø 14 mm
8853142	test tube rack for 42 test tubes Ø 16 mm
8853150	test tube rack for 36 test tubes Ø 18 mm
8853169	test tube rack for 33 test tubes Ø 20 mm
8853185	test tube rack for 18 test tubes Ø 25 mm
8853177	test tube rack for 12 test tubes Ø 30 mm


Racks for centrifuge tubes, with screw cap

8853088	Rack for 42 centrifuge tubes Ø 16 mm
8853096	Rack for 36 centrifuge tubes Ø 18 mm
8853193	Rack for 33 centrifuge tubes Ø 20 mm
8853240	Rack for 12 centrifuge tubes Ø 30 mm


Clamps for Erlenmeyer flasks (for optional mounting onto the universal trays)


8854505	Clamps for Erlenmeyer flasks 25 ml required / possible max. number of clamps for tray type EU - max. 49 pcs. for tray type FU - max. 98 pcs.
8854513	Clamps for Erlenmeyer flasks 50 ml required / possible max. number of clamps for tray type EU - max. 48 pcs. for tray type FU - max. 96 pcs.
8854521	Clamps for Erlenmeyer-flasks 100 ml required / possible max. number of clamps for tray type EU - max. 24 pcs. for tray type FU - max. 48 pcs.
8854556	Clamps for Erlenmeyer flasks 250 ml required / possible max. number of clamps for tray type EU - max. 17 pcs. for tray type FU - max. 39 pcs.

Cat.-no.	Features and Specifications
8854572	Clamps for Erlenmeyer flasks 500 ml required / possible max. number of clamps for tray type EU - max. 12 pcs. for tray type FU - max. 26 pcs.
8854599	Clamps for Erlenmeyer flasks 1.000 ml required / possible max. number of clamps for tray type EU - max. 8 pcs. for tray type FU - max. 17 pcs.
8854610	Clamps for Erlenmeyer flasks 2.000 ml required / possible max. number of clamps for tray type EU - max. 4 pcs. for tray type FU - max. 9 pcs.
8854629	Clamps for Erlenmeyer flasks 3.000 ml, required / possible max. number of clamps for tray type EU - max. 4 pcs. for tray type FU - max. 8 pcs.
8854637	Clamps for Erlenmeyer flasks 5.000 ml required / possible max. number of clamps for tray type EU - max. 2 pcs. for tray type FU - max. 6 pcs.
	Clamps for Fernbach flasks (optional equipment of the universal trays)
8854564	Clamps for Fernbach flasks 450 ml required / possible max. number of clamps for tray type EU - max. 6 pcs. for tray type FU - max. 15 pcs.
8854600	Clamps for Fernbach flasks 1.800 ml required / possible max. number of clamps for tray type EU - max. 1 pc. for tray type FU - max. 6 pcs.
8854640	Clamps for Fernbach flasks 2.800 ml required / possible max. number of clamps for tray type EU - max. 1 pc. for tray type FU - max. 6 pcs.

Cat.-no.	Features and Specifications
	Plastic clamps for Erlenmeyer flasks (optional equipment of the universal trays)
8854700	Plastic clamps for Erlenmeyer flasks 100 ml required / possible max. number of clamps for tray type EU - max. 20 pcs. for tray type FU - max. 58 pcs.
8854711	Plastic clamps for Erlenmeyer flasks 250 ml required / possible max. number of clamps for tray type EU - max. 20 pcs. for tray type FU - max. 40 pcs.
8854722	Plastic clamps for Erlenmeyer flasks 500 ml required / possible max. number of clamps for tray type EU - max. 16 pcs. for tray type FU - max. 26 pcs.
8854733	Plastic clamps for Erlenmeyer flasks 1.000 ml required / possible max. number of clamps for tray type EU - max. 9 pcs. for tray type FU - max. 15 pcs.
Note:	Plastic clamps are ideal for use with Erlenmeyer flasks with tubing connection. See below! Shaking flask GL Erlenmeyer type, with baffles for increased turbulence
Fig. 45: Shaking flask Erlenmeyer type; with baffles in the concave side-wall section and connector GL14 for tubing connection	 <p>shaking flask made of DURAN®-glass, Erlenmeyer-type;; straight neckpiece for metal caps; neckpiece of Ø = 38 mm; flask body with three baffles at ∠ 120deg with GL 14 - connector incl. plastic cap and hose connectors for tubings of size 4 x 7 mm</p>
8861064	size 300 ml, max. Ø = 87 mm, H = 161 mm; set of 10 pcs.
8861072	size 500 ml, max. Ø = 105 mm, H = 183 mm; set of 10 pcs.
8861080	size 1.000 ml, max. Ø = 131 mm, H = 232 mm; set of 10 pcs

Cat.-no.	Features and Specifications
Accessories for baffles shaking flasks	
8861099	cap made of aluminium, set of 10 pcs.
8861102	stainless steel cap, set of 10 pcs.


Cat.-no.	Features and Specifications
Shaking flask Erlenmeyer type, with baffles for increased turbulence	
<p>Fig. 46: Shaking flask Erlenmeyer type; with baffles in the concave side-wall section</p>	 <p>shaking flask made of DURAN[®]-glass, Erlenmeyer-type; straight neckpiece for metal caps; neckpiece of Ø = 38 mm or version with cellucotton stoppers flask body with three baffles at ∠ 120deg</p>
8861005	size 300 ml, max. Ø = 87 mm, H = 161 mm; set of 10 pcs.
8861013	size 500 ml, max. Ø = 105 mm, H = 183 mm; set of 10 pcs.
8861021	size 1.000 ml, max. Ø = 131 mm, H = 232 mm; set of 10 pcs.
8861022	size 2.000 ml, max Ø = 166 mm, H = 305 mm; set of 10 pcs.
8860998	size 500 ml, version for cellucotton stoppers, etc. max Ø = 105 mm, H = 180 mm; set of 10 pcs. Other versions on request

Universal Mounting Set	
Basic element type B,	
each consisting of 2 side frames and 4 rods	
<p>Fig. 47: Basic element type B of universal mounting set</p>	
8854238	type B-2 for trays EU
8854246	type B-3 for tray FU

Cat.-no.	Features and Specifications
----------	-----------------------------


Universal clamping rod type U

Fig. 48:
Universal clamping rod


8854254 type U for basic elements B-2 and B-3

Fig. 49:
Mounting system with
universal clamping rods

Mounting set for separating funnels type S

Fig. 50:
Mounting set for
separating funnels type S

8854262 type S-1 for 5 pcs. separating funnels 50 or 100 ml; set includes clamps,
retaining springs and a universal clamping rod U-2

8854270 type S-2 for 3 pcs. separating funnels 250, 500 or 1.000 ml; set includes clamps,
retaining springs and a universal clamping rod U-2

Max. number of mounting sets for tray type is

1 mounting set S 1 or S 2 for tray type EU with set B-2

2 sets of mountings S 1 or S 2 for tray type FU with set B-3

4.2.4 Additional Laboratory Equipment and Consumbles

Cat.-no.	Ordering data, features and specifications
----------	--

8521620	FRIGOMIX® R circulation cooler cooling capacity 500 W at + 10°C temperature range -10 ... +40°C pump capacity : 16 l/min., 300 mbar version 230 V / 50 Hz
9321802	Meliseptol® alcoholic agent for disinfection and cleaning, ready for use bottle 250 ml, with spray diffuser

- 4.3 Safety Instructions Laboratory Shakers and Incubation Shaking Cabinets**
(Safety Instructions Laboratory Shakers and Incubation Shaking Cabinets see appendix)
- 4.4 EG Declaration of Conformity**
(EG Declaration of Conformity see appendix)
- 4.5 Declaration of Decontamination**
(Declaration of Decontamination see appendix)
- 4.6 Information and Instructions on Disposal and Repairs**
(Information and Instructions on Disposal and Repairs see appendix)
- 4.7 Dimension sheet CERTOMAT® BS-T**
(Dimension sheet see appendix)
- 4.8 Prospectus CERTOMAT® BS-T**
(Prospectus CERTOMAT® BS-T see appendix)

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany

Telephone +49.551.308.0
Telefax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Göttingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission of
Sartorius Stedim Biotech GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology,
features, specifications, and
design of the equipment
without notice.

Status:
March 2008,
Sartorius Stedim Biotech GmbH,
Göttingen, Germany

sartorius stedim
biotech

Safety Instructions

Laboratory Shakers and Incubation Shaking Cabinets

85030-520-67

Content

1	Introduction	3
2	Transport and Installation	3
2.1	Transporting Equipment	3
2.2	Requirements at the site of installation	3
3	Operating Instructions	4
3.1	Working with Media	4
3.2	Initial Startup and Normal Operation	4
4	Cleaning, Maintenance, and Service	5
4.1	Cleaning	5
4.2	Maintenance and Service	5
4.3	Returning Defective Units	5
4.4	Relocating equipment	5

1 Introduction

These Safety Instructions contain our recommendations for precautions in handling CERTOMAT® laboratory shakers and incubation shaking cabinets from the product program of the company Sartorius Stedim Biotech GmbH, especially for

- transport and installation
- equipping and initial startup
- operation
- cleaning, maintenance, and service

These safety instructions are non-transferable and meant exclusively for the equipment delivered. They are not valid for other equipment. Additional, specific safety regulations (i.e., legally or otherwise mandatory rules) are not dealt with by these instructions. Please consult your local authorities for further information.

Possible risks and dangers are marked with this symbol and highlighted like this paragraph. Ignoring these warnings may cause damage to the instrument or lead to other material or personal damage.

Steps to be carried out with special care or special aspects or issues to be considered are marked with safety instructions formatted as this paragraph.

[→ ..] Identifies references to contents of these instructions or other documents. Labels of illustrations, sections, or documents are in parentheses.

2 Transport and Installation

2.1 Transporting Equipment

The shakers and especially the incubation shaking cabinets are heavy. Use only suitable transportation and lifting tools when moving or setting up units.

Make sure clearances are sufficient along the transport route and that the combined weight of lifting and transportation tools does not exceed the maximum load permitted for the floor or ground.

Use only qualified personnel to transport the equipment. Make sure personnel is protected from injury during transport.

Do not remove any transport protection before the equipment has been transported to the site of installation.

2.2 Requirements at the site of installation

1. The floor or the laboratory bench where the equipment is to be installed must be able to support the weight of the fully equipped units and systems.
2. The installation place should be sufficiently dimensioned and feature a non-slip surface.
3. The equipment has to be easily accessible for operation, maintenance, and service tasks.
4. The shakers and the incubation shaking cabinets may transmit vibrations to the installation surface especially when carrying a heavy load. For this reason, install the shakers so that surrounding equipment is not affected.
5. The power supply must match the rating indicated on the type plate and must be grounded correctly.

Use a level to align the equipment after setting it up!

Before stacking equipment (max. 3 units), make sure that the bottom most device is properly aligned horizontally!

3 Operating Instructions

1. Make sure, that only authorized people have access to the working area.
2. Personnel has to be instructed about applicable safety instructions and other legally or otherwise compulsory regulations.

3.1 Working with Media

The utilized media and substances can be hazardous and may pose dangers that are specific to the substances and the processes and thus cannot be described here in detail. Your company should issue corresponding safety regulations and carefully instruct the operating personnel.

Recommended basic protection equipment includes such items as suitable working clothes, gloves, goggles, and respiratory equipment, if necessary.

3.2 Initial Startup and Normal Operation

1. Use only units, equipment, accessories, and spare parts released by the company Sartorius Stedim Biotech GmbH for use with the shaker.
2. Check all components for damage, especially glass components such as jars. Any damaged components must not be used.
3. Never use the shaker without the tray. Switch the unit off before installing/uninstalling trays, mounting sets, or additional accessories.
4. When mounting the trays, leave enough space between the unit and other objects, especially in case of trays extending over or beyond the shaker table or bench. Make sure that trays are properly installed and fastened.
5. Fasten the tray mountings and accessories properly. Distribute loads evenly on the tray when loading the shaker with the shaking vessels (Erlenmeyer flasks, bottles, etc.).
6. In case of incubation shaking cabinets, make sure the tray is firmly positioned; close the doors before starting the shaking process.
7. Before starting the shaker and during normal operation, check that the vessels are firmly positioned in their holders. The vessels must not become loose, rattle, or touch one another or the walls of the shaker. Switch the unit off immediately if vessels become loose in their holders.
8. Operate the equipment only within the permitted range of operating parameters. Comply with the equipment instructions and the technical data.

Head the lateral movement of the trays and the resulting impact or crushing hazards, especially with large trays, in case of low clearances to other equipment or objects or large amplitudes and heavy loads.

Never put fingers between tray support and frame when the unit is active. The space between these components is narrow due to constructional constraints. Risk of injury!

9. After switching off the shaker, wait until the shaker support has come to a complete standstill before removing or adding vessels.

4 Cleaning, Maintenance, and Service

Regular cleaning and maintenance ensures the equipment's functionality and safety.

Comply with the legal safety regulations applicable to the field of application concerning the handling of microbial material and contaminated equipment.

Prior to any cleaning or maintenance work, disinfect or sterilize any biologically contaminated components according to the safety regulations.

Prior to any cleaning or maintenance work, switch the unit off and pull the power cable to disconnect the unit from the mains power supply. Make sure that the unit cannot be restarted or activated accidentally.

Servicing and repairing drive systems, replacing drive belts, or working on electrical/electronic equipment should only be performed by qualified service personnel.

4.1 Cleaning

1. Clean the surfaces of the units with standard laboratory cleaners. Do not use corrosive agents (e.g., chlorous agents).
2. Broken glass or spilled liquids should be removed as quickly as possible.
3. For incubator shakers special hints are given in the operating manual.

4.2 Maintenance and Service

Maintenance and troubleshooting by the user should be restricted to the following tasks:

- Checking the voltage supply in the lab and the equipment connections.
- Checking and replacing defective fuses.
- For incubation shaking cabinets: Replacing defective light fittings.

If problems or malfunctions cannot be solved or remedied, please contact your Sartorius AG representative or contact us directly:

Sartorius AG
Servicezentrum Nord
Weender Landstraße 94-108
37075 Göttingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Fax +49.551.308.3730

A written cleaning and decontamination certificate has to be provided when requesting service technicians of the company Sartorius AG to perform maintenance tasks.

4.3 Returning Defective Units

Defective units or components can be returned to the responsible Sartorius AG representative or directly to Sartorius AG.

1. The equipment has to be clean, hygienically flawless, and carefully packaged.
2. Contaminated parts have to be cleaned, disinfected, or sterilized according to the respective and applicable safety guidelines. Compliance with the guidelines has to be documented with a cleaning and decontamination certificate to be enclosed with the shipment.
3. Negotiate the return shipment with the Service Dept. before shipping any equipment.

4.4 Relocating equipment

1. Comply with the described cleaning and maintenance measures described to avoid biological or chemical hazards.
2. Disconnect the equipment from the mains power supply. Remove vessels and additional equipment or components inside of the unit or connected to the unit.
3. If applicable, reinsert the transport protection. Read and comply with the additional notes under "Transporting Equipment" (2.1).

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany

Telephone +49.551.308.0
Fax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Göttingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission of
Sartorius Stedim Biotech GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology,
features, specifications, and
design of the equipment
without notice.

Status:
February 2008,
Sartorius Stedim Biotech GmbH,
Göttingen, Germany

Printed in Germany on paper that has
been bleached without the use of chlorine
W4A000 · KT
Publication No.: SB6043az08022
Order No.: 85030-520-67

sartorius stedim
biotech

Declaration of Conformity

according to EU Machinery Directive 89/392/EU,

Appendix II A

Company

Sartorius Stedim Biotech GmbH

Address

August-Spindler-Strasse 11, 37079 Goettingen; Germany
Phone +49.551.308.0, Fax +49.551.308.32 89
www.sartorius-stedim.com

We herewith declare that the device described below fulfills the relevant fundamental safety requirements and health regulations specified by the appropriate EU-Directive, with respect to its design and construction and to the version as commercialized.

This declaration becomes legally invalid if modifications are performed on the device which have not been certified by Sartorius Stedim Biotech GmbH.

Designation of the device

CERTOMAT® BS-T

Model, version

Incubation Shaking Cabinet

Cat.-No.

8865426, 8865434, 8865523, 8865531, 8865620
8865639, 8865728, 8865736

Relevant directives of the EC

EU Machinery Directive (89/392/EU) in version 93/44/EU
EU Directive on Electromagnetic Compatibility (89/336/EU) in version 92/31/EU
EU Low-Voltage Equipment Directive (73/23/EU)

Applied harmonized standards

EN 292-1, EN 292-2
EN 50081-1, EN 50082-1, EN 55011
EN 61010-1

Applied national standards and technical specifications

Date and Signature

11.01.08

Function of the Signatory

V. Niebel
Managing Director

Dr. Susanne Gerighausen
Quality Management
Representative

F-13-01 Version 01	Page 1 / 4 Effective Date 01.01.2008	Return of material	
--------------------------	---	---------------------------	--

Declaration about decontamination and cleaning of equipment and components

When returning equipment or components, please describe on page 2 of this form the problem(s) or fault(s) you have found. Please also indicate the remedial actions you require.

To protect our personnel, we require all equipment or components be free of biological, chemical, or radioisotopic contaminants. We will only accept such equipment or components when:

- the equipment or components have been adequately cleaned and decontaminated.
- this declaring document has been completed, signed and returned by an authorized person.

Please help us in assuring a safe, hazard-free work environment.

A. Description of the Equipment or Component(s)

Description / Cat. No.	
Serial no.	
No. of invoice/delivery note	
Date of delivery	

B. Contamination / Cleaning

Attention: Please specify exactly the biological, chemical, or radioisotopic contaminant	Attention: Please describe the cleaning and decontamination procedure/method
The equipment was contaminated with	and it has been cleaned and decontaminated by

F-13-01 Version 01	Page 2 / 4 Effective Date 01.01.2008	Return of material	
------------------------------	---	---------------------------	--

C. Legally binding declaration

I (we) certify that all information given in this form is correct and complete.
The equipment and components have been adequately decontaminated and cleaned according to the legal requirements.
No chemical or biological or radioisotopic risks remain that can endanger exposed persons' safety or health.

Company / Institute	
Address / Country	
Tel. / Fax (with area code)	
Name of the authorized person	
Position	
Signature / Date	

D. Reason for return

☐ wrong delivery
☐ exchange
☐ repair
☐ modification
☐ disposal
☐ other

E. Please describe the problem(s) or fault(s) you have found (for repair) and/or indicate the remedial actions you require

F-13-01 Version 01	Page 3 / 4 Effective Date 01.01.2008	Return of material	 sartorius stedim biotech
--------------------------	---	---------------------------	---

F. Reserved for Sartorius-Service Center

Notes

Please pack the unit properly and send it freight paid to your local service supplier or directly to Sartorius AG Servicezentrum Nord, Germany.

Sartorius AG
Servicezentrum Nord
Weender Landstraße 94-108
37075 Göttingen, Germany
Telephone +49.551.308.3729 / 3740 / 3741
Fax +49.551.308.3730

F-13-01 Version 01	Page 4 / 4 Effective Date 01.01.2008	Return of material	 sartorius stedim biotech
------------------------------	---	---------------------------	---

Sartorius Stedim Biotech GmbH
August-Spindler-Straße 11
37079 Göttingen, Germany

Telephone +49.551.308.0
Telefax +49.551.308.3289
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Göttingen, Germany.
All rights reserved. No part
of this publication may be
reprinted or translated in any
form or by any means without
the prior written permission
of Sartorius Stedim Biotech
GmbH.

The status of the information,
specifications and illustrations in
this manual is indicated by the
date given below. Sartorius Stedim
Biotech GmbH reserves the right
to make changes to the technology,
features, specifications, and
design of the equipment
without notice.


Status:
February 2008,
Sartorius Stedim Biotech GmbH,
Göttingen, Germany

Information and Instructions on Disposal and Repairs

Packaging that is no longer required must be disposed of at the local waste disposal facility. The packaging is made of environmentally friendly materials that can be used as secondary raw materials.

The equipment, including accessories and batteries, does not belong in your regular household waste. The EU legislation requires its Member States to collect electrical and electronic equipment and disposed of it separately from other unsorted municipal waste with the aim of recycling it.

In Germany and many other countries, Sartorius AG takes care of the return and legally compliant disposal of its electrical and electronic equipment on its own. These products may not be placed with the household waste or brought to collection centers run by local public disposal operations – not even by small commercial operators.

For disposal in Germany and in the other Member States of the European Economic Area (EEA), please contact our service technicians on location or our Service Center in Goettingen, Germany:

Sartorius AG
Service Center
Weender Landstrasse 94-108
37075 Goettingen, Germany

In countries that are not members of the European Economic Area (EEA) or where no Sartorius affiliates, subsidiaries, dealers or distributors are located, please contact your local authorities or a commercial disposal operator.


Prior to disposal and/or scrapping of the equipment, any batteries should be removed and disposed of in local collection boxes.

Sartorius AG, its affiliates, subsidiaries, dealers and distributors will not take back equipment contaminated with hazardous materials (ABC contamination) – either for repair or disposal. Please refer to the accompanying leaflet/manual or visit our Internet website (www.sartorius.com) for comprehensive information that includes our service addresses to contact if you plan to send your equipment in for repairs or proper disposal.

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CERTOMAT® BS-T
Programmable
top loading incubation
shaking cabinet

CERTOMAT® BS-T

Programmable top loading incubation shaking cabinet

The CERTOMAT® BS-T is a high volume shaking cabinet for all laboratories in micro-biology, molecular biology as well as bio-chemistry, plant physiology and others.

An innovative programming of all parameters is done via a membrane key pad in real time. This control allows to store up to 5 programs with 5 steps each password protected. Via a repeat function, unlimited repetition is possible.

Process safety is maintained during power failure. Time and duration are stored as an alert. The Memory Function restores the last set points after power failure.

Continuous recording for documentation of all parameters is possible via analog out.

The shaking amplitude of 25 mm and 50 mm can be converted.

A built-in cooling as option allows temperatures near and beyond ambient and is needed for the optional illumination unit.

The interior made of polished stainless steel (1.4301) can easily be kept in hygienic conditions.

The standard spill tray collects liquid from broken glass. Tray fixation is done in one motion, safe and easy.

The innovative programming of the CERTOMAT® BS-T allows individual control of temperature, shaking speed and light.

A pre-step allows a cooling period before the start of the incubation program. Cultures are visible through the large double panned window in the lid.

Temperatures at or down to 10 °C below ambient are reached with an optional integrated cooling unit.

Cultivation of phototrophic organisms is possible with the optional illumination unit installed in the lid. The lights are programmable for simulation of day/night conditions.

Benefits

Programming of all parameters

Stores up to 5 programs

High process safety

Output for continuous recording

Indicates power failure by length and time

Orbit 25 mm and 50 mm can be converted


Stainless steel spill tray to protect from damage

Built-in cooling unit optional

Illumination unit optional

Tray fixation fast and easy



Accepts 6 liter flasks

Zubehör

Reference	Description
886 1463	Illumination unit for CERTOMAT® BS-T, 5×18 W, individually activated, programmable, only in combination with cooling
886 4489	Support frame, welded sectional frame construction, height adjustable feet
Tray Type F (800×420 mm)	
885 3738	74 clamps inox for flasks 100 ml
885 3762	40 clamps inox for flasks 250 ml
885 3789	26 clamps inox for flasks 500 ml
885 3800	15 clamps inox for flasks 1,000 ml
885 3037	Universal tray to be completed with clamps, racks or mounting system Type FU (800×420 mm)
Stainless steel clamps, Capacity for tray type FU	
885 4505	for Erlenmeyer flasks 25 ml (max. 98)
885 4513	for Erlenmeyer flasks 50 ml (max. 96)
885 4521	for Erlenmeyer flasks 100 ml (max. 48)
885 4556	for Erlenmeyer flasks 250 ml (max. 39)
885 4572	for Erlenmeyer flasks 500 ml (max. 26)
885 4599	for Erlenmeyer flasks 1,000 ml (max. 17)
885 4610	for Erlenmeyer flasks 2,000 ml (max. 9)
885 4629	for Erlenmeyer flasks 3,000 ml (max. 8)
885 4637	for Erlenmeyer flasks 5,000 ml (max. 6)
885 4564	for Fernbach of 450 ml (max. 15)
885 4600	for Fernbach of 1,800 ml (max. 6)
885 4640	for Fernbach of 2,800 ml (max. 6)
Plastic clamps reinforced with glass fibre	
885 4700	for Erlenmeyer flasks 100 ml (max. 48)
885 4711	for Erlenmeyer flasks 250 ml (max. 39)
885 4722	for Erlenmeyer flasks 500 ml (max. 26)
885 4733	for Erlenmeyer flasks 1,000 ml (max. 17)
Hinged racks for test tubes (8 racks max.)	
885 3134	for 64 tubes dia. 14 mm
885 3142	for 42 tubes dia. 16 mm
885 3150	for 36 tubes dia. 18 mm
885 3169	for 33 tubes dia. 20 mm
885 3185	for 18 tubes dia. 25 mm
885 3177	for 16 tubes dia. 30 mm
Hinged racks for centrifugation tubes (8 racks max.)	
885 3088	for 42 tubes dia. 16 mm
885 3096	for 36 tubes dia. 18 mm
885 3193	for 33 tubes dia. 20 mm
885 3240	for 16 tubes dia. 30 mm

Reference	Description
885 0321	Holders for microtiter plates, stainless steel for 1 standard 96-well plate or deepwell plate standard plates: max. 12 holders on EU tray, 21 holders on FU tray deepwell plates: max. 9 holders on EU tray, 18 holders on FU tray
886 4497 886 0416 886 4470	Sticky tape for universal trays Standard, roll of 50 m, 30×1 mm Premium, roll of 10 m, 30×1 mm, repeated use Anti-skid layer, 380×450 mm, for individual cut
885 4246 885 4254	Universal mounting system Basic element B-3 for tray FU Clamping rod, type U for mounting systems B-2 und B-3
885 4262 885 4270	Mounting system for separatory funnels (for use with type FU, in combination with type B-3 basic element) Type S-1 for 5 separatory funnels 50 and 100 ml Type S-2 for 3 separatory funnels 250, 500 and 1,000 ml
886 1005 886 1013 886 1021 886 1022	Shaking flasks, DURAN, Erlenmeyer type, 3 baffles at 120°, straight rim Erlenmeyer flasks 300 ml, pack of 10 Erlenmeyer flasks 500 ml, pack of 10 Erlenmeyer flasks 1,000 ml, pack of 10 Erlenmeyer flasks 2,000 ml, pack of 10
886 1064 886 1072 886 1080	Shaking flasks, DURAN, Erlenmeyer type, 3 baffles at 120°, straight rim, connector GL 14 Erlenmeyer flasks 300 ml, pack of 10 Erlenmeyer flasks 500 ml, pack of 10 Erlenmeyer flasks 1,000 ml, pack of 10
886 1099 886 1102	Caps Cap aluminium, pack of 10 Cap inox, pack of 10
886 0998	Shaking flasks, DURAN, Erlenmeyer type, 3 baffles at 120°, narrow neck for plug Erlenmeyer flasks 500 ml, pack of 10

Technical Data

Ordering information

CERTOMAT® BS-T version with circulation/heating (UH)		
230 V/50 Hz	886 5426	CERTOMAT® BS-T/25 mm
230 V/50 Hz	886 5523	CERTOMAT® BS-T/50 mm
115 V/60 Hz	886 5434	CERTOMAT® BS-T/25 mm
115 V/60 Hz	886 5531	CERTOMAT® BS-T/50 mm
CERTOMAT® BS-T version with circulation/heating/cooling (UHK)		
230 V/50 Hz	886 5620	CERTOMAT® BS-T/25 mm
230 V/50 Hz	886 5728	CERTOMAT® BS-T/50 mm
115 V/60 Hz	886 5639	CERTOMAT® BS-T/25 mm
115 V/60 Hz	886 5736	CERTOMAT® BS-T/50 mm

Technical Data

Mechanical	Shaking mode: Amplitude: Drive: Drive control:	orbital 25 mm or 50 mm brushless motor with triple eccentric bearing Microprocessor controlled
Shaking Speed	Range: Adjustment: Indication: Speed accuracy:	40 to 400rpm membrane keypad alphanumeric LCD ±1% of max. value
Temperature	Temperature range: with cooling: Air circulation:	ambient +8 °C to +70 °C ambient -10 °C to +70 °C 180 m³/h
Trays	Recommended trays: Maximum load: Mounting:	Type F/FU (420 × 800 mm) all trays with accessories out of Sartorius program Spring-clip snap mechanism
Programming	Parameters: Timing: Programs: Adjustment: Indication:	Speed, temperature, illumination (optional) date, hour, minute up to 5 programs with cycling Membrane key pad LCD, alphanumeric, 4 lines, 20 characters
Operating Environment	Safety class: Admissible temperature: Relative humidity:	IP 21 from +10 °C to +45 °C max +32 °C with mode cooling from 10% to 90%, not condensing
Documentation	Analog out:	for speed and temperature 0–10 V, optional 0/4–20 mA
Electrical	Power supply: Heating capacity: Cooling capacity:	230 V/50 Hz 650 W 500 W
Dimensions and weight	External (W × H × D): Internal (W × H × D): Weight (without tray):	1150 × 760 × 750 mm 890 × 535 × 595 mm 171 kg

All models are delivered without tray and other accessories

For incubating cells or mixing liquids a tray is needed together with additional accessories to hold shaking flasks, separation funnels or tubes.

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The reproduction, distribution and utilization of this document as well as the communication of its contents to others without explicit authorization is prohibited. Offenders will be held liable for the payment of damages. All rights reserved in the event of the grant of a patent, utility model or design.
Sartorius Stedim Singapore Pte. Ltd.
Sartorius Stedim Australia Pty. Ltd.